

Technical Data Sheet

=====

Item: Alcohol swabs
Item No.: 1131897

1) Material:

- Non-woven material saturated with 70% Isopropyl alcohol

2) Sealing:

Aluminium foil paper pouches

3) LOT-No:

The LOT-No. and expiry date EXP in format DDMMYY are printed on each pouch, on each box and on the outer cartons

- DD indicates the day, for example, 25
- MM indicates the month, for example, 02;
- YY indicates the year, for example, 12;

4) Use:

For skin cleansing as well as for many jobs in the laboratory.

5) Instruction for use:

Rub skin brskly in a circular motion.

CAUTION: Only for external use: in the event of skin irritation or rash, discontinue use. Avoid contact with eyes; if contacted, immediately flood the eye with plenty of water. If ingested, DO NOT induce vomiting, obtain medical attention immediately.

6) Presentation:

folded: appr. 3 x 3 cm

unfolded: appr. 6 x 3 cm

7) Quantity:

- swabs per box:

100 pcs.

- swabs per carton:

2000 pcs.



© www.herenz.de

www.herenz.de

Phone +49-040-739-20-40
Fax +49-040-730 41 48
Email hhh@herenz.de
ILN 40 49434 00000 8
Alkoholtupfer-e.docx

Managing Directors

Dieter Cyll, Antje Herenz-Cyll
Comm. Reg. Hamburg HRB 48884
Tax No. 44/729/00248
V.A.T. No. DE 118 665 646

Postbank Hamburg

Sort Code 200 100 20
Account No. 6846 46 200
IBAN DE97200100200684646200
Swift-Code PBNKDEFF

Deutsche Bank AG

Sort Code 200 700 00
Account No. 55 00 251
IBAN DE82200700000550025100
Swift-Code DEUTDEHH

Hamburger Sparkasse

Sort Code 200 505 50
Account No. 1085 243 002
IBAN DE87200505501085243002
Swift-Code HASPDEHH